

36th South African
NATIONAL BLOOD
Transfusion Congress

22 – 25 August 2022 - Durban

SHAPING A SUSTAINABLE FUTURE

Establishment of a stable microbiology proficiency testing matrix in transfusion microbiology

Presenter:

Xoliswa Mpumlwana¹, Ute Jentsch¹, Winnie Kruger²

¹SANBS

²University of Johannesburg



Introduction

- Quality control department in SANBS screens blood products in microbiology using bacterial culture methodologies.
- There is no formal PT programme available for these methods.
- Aim of the study
 - ❖ To develop and evaluate a stable matrix-equivalent microbiology PT programme using pooled platelets (PP) and red blood cells (RBC) as the basis matrix.



BacT/ALERT 3D (Culture)



BD FX400 (Culture)




Vitek II compact (ID and AST)

Study design and methodology


- Approval was obtained from SANBS and University of Johannesburg.
- A prospective cross-sectional study approached in 2 phases;
 - ❖ matrix preparation and preparation of samples for testing
 - ❖ sample testing .
- Samples were tested in 2 stages over 3 months.
 - ❖ Testing by the researcher to establish reference results and stability of the PT matrix.
 - ❖ Testing by the participating laboratories for evaluation purposes.
- Methods: Gram staining, Bacterial culture, Bacterial identification and Antimicrobial Susceptibility Testing (AST) in both testing stages.
- Bacterial culture was performed using the BacT/ALERT 3D and Bactec 9120 while Vitek II compact was used for AST.
- Cohen's kappa was used for statistical analysis (>0.75 indicating the high level of agreement).



Study design and methodology



Study design and methodology



RR - Reference results



ST- Stability testing

PL - Participating laboratories



Results by researcher

Stability testing



- A high level of agreement for bacterial culture and Gram stain results was obtained throughout all stability detection days with kappa scores of >0.75
- Day 1 = Reference results



Results by researcher

Stability testing

Bacterial identification stability testing	
Testing days	Kappa value
Day 1 versus Day 3	1,000
Day 1 versus Day 5	1,000
Day 1 versus Day 9	1,000
Day 1 versus Day 14	0,955
Day 1 versus Day 23	0,955
Day 1 versus Day 30	0,842

- A kappa score of 1.000 was obtained on AST results for all bacteria that required AST testing.



Results by participating laboratories

- Kappa score of 1.000 was obtained on AST results for all bacteria that required AST.

Participants	Bacterial Culture	Gram Stain	Bacterial Identification
Lab 1 vs Reference results	1.000	0.947	0.958
Lab 2 vs Reference results	1.000	0.920	0.958
Lab 3 vs Reference results	1.000	0.920	0.894

- Bacterial culture had 156/156 = 100% compliance
- Gram stain had 151/156 = 97% compliance
- Bacterial Identification had 149/156 = 96% compliance
 - ❖ *Staphylococcus epidermidis*
 - ❖ *Bacillus cereus*
 - ❖ *Staphylococcus aureus*
 - ❖ *Listeria monocytogenes*
 - ❖ *Streptococcus pyogenes*
- AST had 156/156 = 100%



Discussion

Stability results

- PT matrices were stable for up to 23 days in all storage conditions between 2°C and 37°C.
- The results obtained during stability testing demonstrated that storage time of up to 23 days was acceptable, except for *Listeria monocytogenes*.
- *Listeria monocytogenes* was stable for 14 days.
- This strains requires certain growth substances which were not provided in the platelet and RBC concentrates.
- The growth substances include amino acids such as leucine and vitamins such as riboflavin.
- *Morganella morganii* could not be identified on day 30.
- This strain was validated by the Paul-Ehrlich-Institut (PEI) for only 7 days in platelets and never validated in red blood cells.



Discussion

Evaluation from participating laboratories

- All positive and negative cases were identified correctly by the participating laboratories on bacterial culture, showing a high level of agreement.
- There were no discrepancies observed between the reference results and the participating laboratories for AST results, 100% compliance was obtained.
- 3% percent of discrepant results were obtained for Gram stain and bacterial identification
- *Listeria monocytogenes* could not be identified by all 3 participants.



Conclusion

- Blood product based matrices provided a possible solution in closing the gap of unavailable PT matrix for transfusion microbiology
- PT matrices are suitable to be sent to areas where delivery can be achieved within seven days from dispatch.
- This matrix covers Gram staining, bacterial screening, bacterial identification as well as antimicrobial susceptibility testing.
- The matrices need to be replicated and validated to see if suitable for shipments to African Blood transfusion services to support blood safety.



Acknowledgements

Many thanks to the following for all the participation, support and encouragement:

- Winnie Kruger from University of Johannesburg,
- Ute Jentsch from the South African National Blood Service (SANBS).
- National Health Laboratory Service (NHLS) infection control.
- SANBS Quality Control (QC) management.
- SANBS Quality Control and Proficiency Testing (PT) teams.
- Research completed as part of my Master's dissertation.
- Contact information:
 - ❖ QualityControl@sanbs.org.za
 - ❖ +27 (11) 761 9220/ 9225



THANK YOU